




StealthStation® i7™ Integrated Navigation System

Optimize Your Surgical Space



- Versatile integrated navigation
- Flexible equipment interfaces
- Professional services



StealthStation® i7™ Integrated Navigation System



www.medtronicnavigation.com

Medtronic Navigation
826 Coal Creek Circle
Louisville, CO 80027
USA
Toll Free: 877-242-9504
Telephone: 720-890-3200
Fax: 720-890-3500

© 2009 Medtronic Navigation, Inc.
All Rights reserved
Printed in the USA

The following are registered or unregistered trademarks of Medtronic, Inc. and its affiliated companies: AxiEM™, iOR®, O-arm®, PoleStar®, StealthStation® i7™, StealthViz™, Synergy®



Versatile Integrated Navigation

- Minimize equipment congestion, to improve procedural workflow
- Ensure clear and quick access to navigation information through high-definition, wide-screen monitors with touch control capabilities
- Incorporate intra-operative imaging information for confidence before you close
- Incorporate advanced visualization and planning capabilities with StealthViz™ application with StealthDTI™ module
- Adapt to your workflow with tailored Synergy® navigation software

Flexible Equipment Interfaces

- Integrate easily with booms, audio/visual, intra-operative imaging, and other OR components
- Work with the technology you select from your preferred partners or choose an iOR® Suite by Medtronic

Professional Services

- Support you from initial planning to final installation to create the ideal OR solution
- Consult with you to analyze surgical workflow, patient positioning, OR space management, ergonomics, and clinical specialty needs
- Collaborate with architects, hospital staff, and other partners to ensure on-time deployment and verification of your final solution
- Offer comprehensive training, ongoing support, and case coverage to ensure your system runs smoothly

System Specifications

- Intra-operative imaging interfaces – PoleStar® Surgical MRI System, iMRI, iCT, O-arm® imaging system, C-arm
- OR suite interfaces – booms, lights, audio, video, endoscopy, microscopes, PACS, cameras, patient vitals
- User interfaces – navigation integration to central equipment control with Audio/Visual vendors
- Professional services – resident outreach and training, accredited courses, on-site training, customized clinical and technical support, and preventive maintenance
- Hi-definition wide-screen monitors:
 - 61 cm (24") or 48cm (19") Nurses Station Monitors
 - 61 cm (24") Boom-Mounted Monitors (touch and non-touch)
 - Resolution of 1920 x 1200, 60 Hz
- AxiEM™ system compatible – Industry's leading electromagnetic option for neuronavigation
- Microscope interface integration with leading manufacturers
- Transfer of exams to StealthStation® i7™ system via DICOM Query/Retrieve, DICOM Receive, USB stick, or CD/DVD
- Improved camera designed with a greater field of view, a longer camera handle and a laser button for easier aiming and positioning